Tokai Microcatheter

Carnelian

* B M P 0 8 7 *

Carnelian
Support

CE 1639

INSTRUCTIONS FOR USE

English.....	2	Lietuviškai.....	77
Български.....	7	Norsk.....	82
Hrvatski.....	13	Polski.....	87
Česky.....	18	Português.....	93
Dansk.....	23	Română.....	98
Nederlands.....	28	Русский.....	103
Eesti.....	34	Srpski.....	110
Suomi.....	39	Slovensky.....	115
Français.....	44	Slovenščina.....	120
Deutsch.....	50	Español.....	125
Ελληνικά.....	56	Svenska.....	131
Magyar.....	62	Türkçe.....	136
Italiano.....	67	Information on packaging.....	141
Latviski.....	72		

TOKAI MEDICAL PRODUCTS, INC.
1485 Sarayashiki, Taraga-cho, Kasugai-city, Aichi
486-0808 Japan
Tel: +81-568-81-7954
Fax: +81-568-81-7785

Emergo Europe
Westervoortsedijk 60,
6827 AT Arnhem
The Netherlands

TOKAI MEDICAL PRODUCTS

2025/07(Rev.2)

ENGLISH TOKAI MICROCATHETER / INSTRUCTIONS FOR USE

DESCRIPTION

The Tokai Microcatheter “Carnelian” and “Carnelian Support” are catheters with a single lumen capable of accommodating guidewires which are smaller than the indicated size. A hub at the proximal end of the catheter incorporates a strain relief and a luer-lock adapter to facilitate the attachment of the accessories. A radiopaque marker at the distal tip facilitates fluoroscopic visualization. The outer diameter of the catheter is coated with a hydrophilic surface to enhance deliverability to the target vasculature.

Model No.	Catalog No.	Catheter Tip Shape	Outer Diameter Proximal/ Distal Fr (mm)	Inner Diameter Proximal/ Distal Inch (mm)	Total length (mm)	Usable Length (cm)	Hydrophilic Coating Length (mm)
Carnelian SI *Two-marker type	MCC110S-E	Straight	2.8/1.9 (0.93/0.65)	0.021/0.021 (0.52/0.53)	1180	110	600
	MCC125S-E	Straight			1330	125	600
	MCC135S-E	Straight			1430	135	600
	MCC110M-E	Mini 45° angle			1180	110	600
	MCC135M-E	Mini 45° angle			1430	135	600
	MCC110W-E	Double angle			1180	110	600
	MCC135W-E	Double angle			1430	135	600
	MCCW130-E*	Straight			1380	130	600
	MCCW150-E*	Straight			1580	150	1100
	MCC150-E	Straight			1580	150	1100
Carnelian 1.7	MSVS110-E	Straight	2.8/1.7 (0.95/0.57)	0.023/0.017 (0.58/0.43)	1180	110	600
	MSVS110M-E	Mini 45° angle			1180	110	
	MSVS110W-E	Double angle			1180	110	
	MSVS135-E	Straight			1430	135	
	MSVS135M-E	Mini 45° angle			1430	135	
	MSVS135W-E	Double angle			1430	135	
Carnelian 1.9	MSB110-E	Straight	2.8/1.9 (0.94/0.63)	0.026/0.020 (0.65/0.50)	1180	110	600
	MSB110M-E	Mini 45° angle			1180	110	
	MSB110W-E	Double angle			1180	110	
	MSB135-E	Straight			1430	135	
	MSB135M-E	Mini 45° angle			1430	135	
	MSB135W-E	Double angle			1430	135	
Carnelian NT	MXN135X-E	Straight	1.9/1.9 (0.63/0.63)	0.017/0.017 (0.43/0.43)	1430	135	800
	MXN155X-E				1630	155	1100
	MXNS135X-E		1.8/1.6 (0.62/0.56)	0.0155/0.0155 (0.40/0.40)	1430	135	800
	MXNS155X-E				1630	155	1100

Model No.	Catalog No.	Catheter Tip Shape	Outer Diameter Proximal/ Distal Fr (mm)	Inner Diameter Proximal/ Distal Inch (mm)	Total length (mm)	Usable Length (cm)	Hydrophilic Coating Length (mm)
Carnelian 2.2	MCV110-E	Straight	2.9/2.2 (0.95/0.75)	0.027/0.020 (0.69/0.51)	1180	110	800
	MCV110S-E	Straight			1180	110	600
	MCV110M-E	45° angle			1180	110	600
	MCV110Q-E	Cobra angle 1			1180	110	800
	MCV135-E	Straight			1430	135	800
	MCV135S-E	Straight			1430	135	600
	MCV135M-E	45° angle			1430	135	600
	MCV135Q-E	Cobra angle 1			1430	135	800
	MCV150-E	Straight			1580	150	1100
	MCV150M-E	45° angle			1580	150	1100
Carnelian 1.8	MSV110-E	Straight	2.7/1.8 (0.90/0.61)	0.023/0.017 (0.58/0.43)	1180	110	800
	MSV110S-E	Straight			1180	110	600
	MSV110M-E	45° angle			1180	110	600
	MSV110Q-E	Cobra angle 1			1180	110	800
	MSV135-E	Straight			1430	135	800
	MSV135S-E	Straight			1430	135	600
	MSV135M-E	45° angle			1430	135	600
	MSV135Q-E	Cobra angle 1			1430	135	800
	MSV150-E	Straight			1580	150	1100
	MSV150M-E	45° angle			1580	150	1100
	MSV135X-E	Straight			1430	135	600
Carnelian 2.7	MCH105-E	Straight	2.9/2.7 (0.95/0.93)	0.027/0.027 (0.69/0.69)	1130	105	600
	MCH125-E	Straight			1330	125	
Carnelian Support 14	MCP090A-E	Straight	2.6/1.8 (0.85/0.61)	0.021/0.017 (0.54/0.43)	980	90	600
	MCP110A-E	Straight	2.6/1.8 (0.87/0.61)		1180	110	800
	MCP135A-E	Straight			1430	135	1000
	MCP150A-E	Straight			1580	150	1200
Carnelian Support HD	MCPH090A-E	Straight	2.6/1.8 (0.87/0.61)	0.0205/0.017 (0.52/0.43)	980	90	600
	MCPH135A-E	Straight			1430	135	1000
	MCPH150A-E	Straight			1580	150	1200
Carnelian Support DP	MCPMS60A-E	Straight	2.6/1.8 (0.87/0.61)	0.0205/0.017 (0.52/0.43)	680	60	400
	MCPMS80A-E	Straight			880	80	600
Carnelian Support BTA	MCPB150-E	Straight	2.6/1.6 (0.87/0.56)	0.021/0.0155 (0.54/0.40)	1580	150	1200
Carnelian Support 18	MEP070-E	Straight	2.6/1.8 (0.85/0.61)	0.021/0.0195 (0.54/0.50)	780	70	400
	MEP090-E	Straight			980	90	600
	MEP135-E	Straight	2.6/1.8 (0.87/0.61)		1430	135	1000
	MEP150-E	Straight			1580	150	1200

INTENDED USE

The Tokai Microcatheter is intended to be used for infusion of diagnostic, embolic, or therapeutic agents into peripheral vasculature (excluding central circulatory system and neurovasculature), and exchange/support of guidewires.

In addition, some device references are used to guide an intravascular medical device to a desired target vessel (Applicable to Catalog No. MCH105-E, MCH125-E).

INDICATIONS

The Tokai Microcatheter is indicated to treat patients with any disease requiring percutaneous catheter introduction for the treatment of the following:

- Diseases that require infusion of diagnostic, embolic, or therapeutic agents
- Diseases that require exchange/support of guide wires

CONTRAINDICATIONS

- Contraindicated in cases where such a procedure could compromise patient health, according to the medical judgment of a physician.
- Contraindicated for use in newborns, pediatric patients, and pregnant women.

TARGET TREATED POPULATION

The Tokai Microcatheter is intended for use in adults, regardless of sex, who are not part of a contraindicated group.

THE INTENDED USERS

This catheter should be used only by physicians thoroughly trained and experienced in percutaneous, intravascular techniques and procedures.

CLINICAL BENEFITS

The Tokai Microcatheter does not have a direct therapeutic or diagnostic function itself. Therefore, its clinical benefits are indirect. Due to device performance characteristics and indications described by TOKAI, the Tokai Microcatheter allows for the following:

- Successful injection of diagnostic/embolic/therapeutic agents
- Successful exchange of guide wire/ Successful passage of guide wire through lesions (Carnelian Support series)

WARNINGS

- Contents are supplied sterile
Do not use if sterile packaging is unintentionally opened or damaged. Inspect the package for any fault or damage to ensure the sterility of the device.
- Do not use the microcatheter in case of any surface irregularities, bends, or kinks. Any damage to the microcatheter may affect its performance and change its characteristics.
- Use the catheter before the Use By date indicated on the package.
- Single-use only
Reuse, re-processing, or re-sterilizing may also create a risk of contamination of the device, cause infection or cross infection, or may be the source of an infectious disease being transmitted from one patient to another. Contamination of the device may lead to injury, illness, or death.
- Exercise care in handling of the catheter during a procedure. If any resistance is met, discontinue the procedure and identify the cause of the resistance by fluoroscopy. Movement of the microcatheter or guidewire against resistance may result in separation of the microcatheter or guidewire tip, damage to the microcatheter or guidewire tip, or vessel perforation.
- The infusion pressure with this catheter should not exceed the maximum pressure. Infusion pressure in excess of the maximum pressure may result in catheter rupture or device failure.
Furthermore, Infusion pressure in excess of the maximum pressure may result in patients or users injury. If the flow of material through the catheter becomes restricted, do not attempt to clear the catheter lumen by infusion. Identify the cause of the blockage and take corrective action, or replace the catheter with a new one before resuming infusion.
- When infusing a drug or contrast medium, the catheter tip should be monitored under fluoroscopy. If no flow is observed, discontinue the infusion of materials, and replace the catheter with a new

one. When the catheter lumen is obstructed, it may cause catheter expansion, rupture, or cracking resulting in leakage of materials even if the infusion pressure is lower than the maximum pressure.

- Comply with the instruction manuals for medical devices and drugs used together with the catheter.
- Before using the catheter, flush the lumen and outer surface with sterile heparinized saline to activate the hydrophilic coating. It is recommended that flush the lumen with sterile heparinized saline as needed to maintain lubrication during the procedure. If the catheter is used while the outer surface is poorly lubricated, the hydrophilic polymer coating will not function properly to maintain lubrication and the coating may come off during use.
- Overtightening of a hemostasis valve may cause damage to the catheter.
- Do not use in case of lesions with severe calcification.
- When operating the catheter, always ensure that the guide wire is ahead of the catheter.
- Do not use the Tokai Microcatheter for any longer than 24 hours.

PRECAUTIONS

- If this catheter is inserted in a catheter with a stopcock, operating the stopcock may cause catheter rupture or device failure.
- Immersing this catheter in or wiping this catheter with chemicals containing organic solvents such as alcohol may cause damage to the catheter, or may compromise lubrication.
- Using the attached syringe for injecting oil-based contrast medium, fat emulsion, or preparations containing fat emulsion may cause damage to parts made of polycarbonate.
- When the catheter is used to replace or support the guidewire, do not rotate the tip of the catheter to prevent damage or disconnection.
- When using coiled embolization materials, select proper ones by checking the sizes of this product and aids in order to ensure that they fit together correctly.
- If embolic materials, etc. stick to the catheter, discontinue use and replace with a new catheter.
- When removing the catheter from the hoop, remove the catheter carefully to avoid kinking.

DISPOSAL

There is a risk of infection if you come into direct contact with a device that has come into contact with blood.

Care should be taken when handling the tip of both the guide wires and the mandrel used with this device, as they can break through the protective gloves

Wind the device properly, being careful not to splash blood. Place in a disposal bag that does not allow the device to stick out. Ask a company specializing in industrial waste to dispose of it properly in accordance with the laws and regulations of each country.

<Precautions for Use>

- 1) Thoroughly flush the catheter with sterile heparinized saline solution prior to use.
- 2) Use the catheter immediately after the sterile packaging is opened. Discard the product after use.
- 3) In case of a shaped product, be careful not to pinch the shaft with your fingers or the like strongly when pulling out the mandrel from the product. Also, if it is hard to pull out the mandrel, do not try to pull it out forcibly but gently little by little.
- 4) Handle with care to prevent accidental damage, bending, or kinking of the catheter when it is inserted into the guide catheter.

COMPLICATIONS

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Possible complications include, but are not limited to, the following:

- Infection
- Angiospasm
- Hepatic failure
- Renal failure
- Ascites/pleural effusion
- Non-target embolism
- Abscess
- Hematoma

- Damage to vessel wall
- Artery dissection /Aneurysm
- Perforation of vessel or arterial wall
- Death
- Post-embolization syndrome
- Pulmonary embolism due to coil migration
- Thromboembolism
- Allergic reaction
- Bleeding

INSTRUCTIONS FOR USE

- 1) Remove the microcatheter with hoop from the sterile packaging.
- 2) Flush the hoop with heparinized saline to activate the hydrophilic coating. The luer fitting attached to the hoop may facilitate the flushing of the hoop.
- 3) Remove the microcatheter from the hoop, and confirm that the surface is well lubricated. If the microcatheter is difficult to remove, repeat injection or place in heparinized saline bath.
- 4) After removing the microcatheter from the hoop, thoroughly flush the lumen with sterile heparinized saline by attaching saline-filled syringe to the catheter hub. Use the attached syringe only for the purpose of cleaning the catheter lumen; do not use it for infusing drug or contrast medium into the body. Do not overtighten the syringe when connecting it to the hub.
- 5) Carefully Insert and advance an appropriate guidewire into the microcatheter.
- 6) Pass the microcatheter wire assembly unit through the hemostasis valve and carefully insert into the appropriate guide catheter. Verify the product size and select an appropriately sized guide catheter. Carefully tighten the valve around the microcatheter to prevent backflow, allowing some movement through the valve by the microcatheter.
- 7) Advance the guidewire and the microcatheter to a selected vascular site by alternately advancing the guidewire and then tracking the microcatheter over the guidewire under fluoroscopy.
- 8) When using a triaxial system, carefully install the insertable microcatheter wire assembly unit before inserting the microcatheter (if applicable) into the guide catheter, or after inserting the microcatheter into the target site.
Note) Compatible sizes: External microcatheter lumen 0.027 inch (0.69 mm) or larger, internal microcatheter outer diameter 1.9 Fr (0.63 mm) or smaller
- 9) If reinserting the guide wire or microcatheter, check to ensure there are no surface irregularities, bends, or kinks. During treatment, flush the lumen with sterilized heparinized saline as needed to maintain lubrication.
- 10) To infuse, completely remove the guidewire from the microcatheter. Connect a syringe with infusate or introduce embolic materials to the microcatheter hub and infuse as required.
- 11) Wind the device properly, being careful not to splash blood. Place in a disposal bag that does not allow the device to stick out. Ask a company specializing in industrial waste to dispose of it properly in accordance with the laws and regulations of each country.

Tokai Medical Products, Inc. assumes no liability with respect to devices reused, re-processed, or resterilized.

Tokai Medical Products, Inc. makes no warranties, expressed or implied, including, but not limited to, merchantability and fitness to the intended use with respect to such device.

Any serious incident associated with use of this medical device should be reported to the device manufacturer and Competent Authority of the country where it occurred.

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